Phenotiki – True phenotyping-in-a-box solution User's manual

Massimo Minervini¹, Mario Valerio Giuffrida¹, Sotirios A. Tsaftaris² ¹IMT Institute for Advanced Studies Lucca, Italy ²University of Edinburgh, UK Manual Version 0.2

November 21, 2016

Contents

1	Inst	allation	5
	1.1	Run source code	5
	1.2	Install Phenotiki from binaries	6
2	Get	ting started with Phenotiki	7
	2.1	Pot Tray Analysis	7
	2.2	Leaf Annotation	13
	2.3	Leaf Counting	15
	2.4	Data Extraction	18
3	Tro	ubleshooting	20
	3.1	General	20
	3.2	Operating system-dependent issues	20
	3.3	Pot Tray Analysis	21
	3.4	Leaf Annotation	22
	3.5	Leaf Counting	22
	3.6	Phenotiki License	22

Disclaimer



Software of the *Phenotiki* project is provided 'AS IS' without warranty of any kind, either express or implied, including, but not limited to, the implied warranties of fitness for a purpose, or the warranty of non-infringement. Without limiting the foregoing, the authors make no warranty that:

- i. the software will meet your requirements
- ii. the software will be uninterrupted, timely, secure or error-free
- iii. the results that may be obtained from the use of the software will be effective, accurate or reliable
- iv. the quality of the software will meet your expectations
- v. any errors in the software obtained from the authors will be corrected.

Software and its documentation made available on the Phenotiki website:

- i. could include technical or other mistakes, inaccuracies or typographical errors. The authors may make changes to the software or documentation made available on its website.
- ii. may be out of date, and the authors make no commitment to update such materials.

The authors assume no responsibility for errors or omissions in the software or documentation available from its website.

In no event shall the authors be liable to you or any third parties for any special, punitive, incidental, indirect or consequential damages of any kind, or any damages whatsoever, including, without limitation, those resulting from loss of use, data or profits, whether or not the authors have been advised of the possibility of such damages, and on any theory of liability, arising out of or in connection with the use of this software.

The use of the software obtained by the authors is done at your own discretion and risk and with agreement that you will be solely responsible for any damage to your computer system or loss of data that results from such activities. No advice or information, whether oral or written, obtained by you from the authors or from the Phenotiki website shall create any warranty for the software.

Changelog

Version 0.2: 05/11/2016

- Pot Tray Analysis:
 - Updated Section 2.1, paragraph Load image from Phenotiki Device, adding the new import methods implemented in the software;
 - Added the paragraph "Extract Masks";
- Leaf Counting:
 - Added Figure 2 and relative text;
 - Updated Figure 6 and 8 and relative text;
 - Updated Table 4 with new SVR, RF, and Time Features paramters;
- Troubleshooting:
 - Added Section 3.2 Operating system-dependent issues;
 - Added a new Q&A on some problems users experienced trying to install *libsvm*.

Preamble

This manual explains the usage of the *Phenotiki* software. It will be assumed that you have acquired images using the Phenotiki acquisition device (further details on www.phenotiki.com). Nonetheless, this software allows you to analyse time-lapse images of top-view rosette plants obtained with a different acquisition pipeline. Our software is written in MATLAB and we will assume as reference version R2015b.

We are making every effort to fully document our software. However, the current version of the manual has an entry-level guide to the software. Extensive usage instructions are provided in this manual and contextual help can be obtained within the software by pressing on this icon:



1 Installation

Our software is released as stand-alone application pre-packaged for several platforms and as source code (refer to www.phenotiki.com for further details).

1.1 Run source code

Once you have unpacked the archive file containing the source code of Phenotiki, from MATLAB:

- navigate to the folder containing the unpacked source code (e. g., ~/Desktop/Phenotiki);
- 2. run the file **PhenotikiMain.m**.
- $! \rightarrow$ You might be asked to install the following dependencies: vl_feat [1] and libsvm [2].¹

In MATLAB, navigate to the vl_feat folder, go to the *toolbox* folder, and run the file **vl_setup.m**. To check if the installation was successful, navigate to a different folder and type in the command window **vl_demo**.

 $! \rightarrow$ Installing vl_feat will add entries to your MATLAB search paths. Therefore, we recommend to save the new pathdef.m, in order to avoid to repeat the installation of vl_feat every time you run the Phenotiki software. You can do it by typing the **pathtool** command in your command window and pressing the "*Save*" button.²

Install libsvm: Go to https://www.csie.ntu.edu.tw/~cjlin/libsvm/ and follow the instructions to download the libsvm source code.

- If you have Windows (any supported version by MATLAB) 64-bit:
 - 1. type **pathtool**;
 - 2. press "Add Folder..." button;
 - 3. navigate inside the libsvm folder and select the **windows** subfolder;
 - 4. press "Save" button.
- For other systems (Windows 32-bit, Linux 32/64-bit, Mac OS X)³:
 - 1. in MATLAB, navigate to the libsvm folder, then go to *matlab* sub-folder;
 - 2. type make;
 - 3. if no errors occurred, type **pathtool**;
 - 4. press "Add Folder..." button;
 - 5. navigate to the libsvm folder and select windows sub-folder;
 - 6. press "Save" button.

¹ Suggested versions: vl_feat 0.9.20 and libsvm 3.20

² On certain systems (e.g., Linux, Mac OS X), it might be required to have administrator (i.e. superuser) privileges to overwrite the main pathdef.m file.

³Learn more about supported compilers of MATLAB R2015b: http:// www.mathworks.com/support/sysreq/files/SystemRequirements-Release2015b_SupportedCompilers.pdf

1.2 Install Phenotiki from binaries

Stand-alone distributions of the Phenotiki software are available for the following operating systems: Windows, Linux, and Mac OS X. The standalone version of our software requires the **MATLAB Compiler Runtime**⁴ (MCR) installed on your system.

 $! \rightarrow$ We kindly recommend to install MCR 2015b.

In case the MCR is not already on your system, the Phenotiki installer process will download and install it. For this step, you are required to have an active Internet connection during the installation.

- $! \rightarrow$ Superuser permissions may be required to complete this step.
 - 1. Launch the installation executable file
 - On Windows: Phenotiki.exe
 - On Mac OS X: Phenotiki.app
 - On Linux: using the terminal, launch the command \$ sudo ./Phenotiki.install
 - 2. in the installation program press "*Next* >";
 - choose the installation directory (leave it as is if you are not sure) and press "Next >";
 - 4. choose the installation directory for the MCR (strongly advised to leave it as is) and press "*Next* >";
 - 5. accept the License Agreement and press "Next >";
 - 6. a summary window will appear, press "Install >";
 - 7. when the installation is complete, press "Finish".

You will find Phenotiki inside the folder specified at step 2.

1.2.1 Set environment variables (only for Linux users)

When the installer finishes, it is required to set two environment variables. Keeping in mind that MCR is installed by default (step 4) in the directory /usr/local/MATLAB/MATLAB_Runtime, add the following lines at the very end of the file /etc/bash.bashrc⁵:

export XAPPLRESDIR=/usr/local/MATLAB/MATLAB_Runtime/v90/X11/appdefaults

export LD_LIBRARY_PATH="\${LD_LIBRARY_PATH}:/usr/local/MATLAB/ MATLAB_Runtime/v90/runtime/glnxa64:/usr/local/MATLAB/ MATLAB_Runtime/v90/bin/glnxa64:/usr/local/MATLAB/ MATLAB_Runtime/v90/sys/os/glnxa64:/usr/local/MATLAB/ MATLAB_Runtime/v90/sys/opengl/lib/glnxa64"

Once you have done it, log-out and log-in back. The software is installed by default in the directory /usr/local/Phenotiki. Run the command \$ /usr/local/Phenotiki/application/Phenotiki to launch Phenotiki.

⁴ Further information at http://www.mathworks.com/products/compiler/mcr/ index.html

⁵ This is valid for Debian and Ubuntu. For different distribution, please check the correspondent file online.



Figure 1: Main window of the Phenotiki Analysis software.

2 Getting started with Phenotiki

Once you have obtained and installed a copy of the Phenotiki software, you are ready to launch it. A start-up window will appear, as displayed in Figure 1. From here, you will have access to all available (currently four) analysis modules:

- **Pot Tray Analysis:** this module will detect and delineate plants in top-view image sequences, separating them from the background (e. g., soil, pot, moss);
- Leaf Labeling: Phenotiki incorporates a leaf annotation tool [3] (available also as a separate software tool at www.phenotiki.com), to delineate each individual leaf of a rosette plant in a computer-aided fashion;
- Leaf Counting: machine learning-based algorithm for counting the number of leaves of rosette plants;
- **Data Extraction:** this module extracts quantitative plant trait descriptors based on the vision tasks performed by the other modules. Extracted data can be imported by external software (e.g., Matlab or R) for further statistical analysis.

Each module is an independent piece of software and you are allowed to run one of them at the time. When you have finished with a module, you need to close the window to return to the main menu.

2.1 Pot Tray Analysis

This module assumes as input a sequence of images showing a top view on a population of rosette plants. Execution of this module will produce the following results:

1. detection and segmentation (i.e. delineation) of each individual plant (cf. Figure 3);



Figure 2: Pot Tray Analysis Tool.

2. extraction of the phenotyping traits listed in Table 1.

We recommend the following workflow: *image loading, plant centres positioning, parameter settings, extracting masks ,* and *getting the traits.*

- **Interface:** Figure 2 shows the interface of the Pot Tray Analysis tool. Specifically, the window is arranged in four panels:
 - A. **File list:** once you load you data set (please see below to see how to do it), this panel will list all the image files with the time stamp of the picture. Through the first column, you can tick which image you want to analyse. This will allow you to run fast experiments on smaller subset of images.
 - B. **Display:** this is the main view point of the interface, where you can see the current time lapse image. You may change it through the horizontal scroll bar or by clicking on the wanted row in the Panel A.
 - C. **Plant list:** once you set the center of the plants (please see below how to do it), you will find the list of the plants with the [x,y] coordinates. Furthermore, if plants belong to different group/genotype/treatment, you may specify it in the second column of this list. At the end of this panel, you find a menu, where you can set the visualization type:
 - *Raw Image:* the original RGB image is displayed;
 - *Detected Plants:* detected plants are showed with different colors (cf. right side of Figure 3);
 - *FG Mask:* Similar as above, but plants are all white;
 - *Contour:* Contour marking the plants are showed on the RGB image.
 - D. *Toolbox:* From this panel you will access to the main functionality of the Pot Tray Analysis. There are four buttons

- i. *Settings:* you can change settings for the segmentation tool (cf. Figure 5);.
- ii. *Extract Mask:* Extract the segmentation mask from the selected plants (cf. Figure 2A).
- iii. *Get Traits:* Extract morphological plant phenotyping traits (cf. Table 1).
- iv. *Save:* Save your progress in a dataset. This file contains all the data of your dataset of images and, once you have it, can be loaded using *Load Dataset* button (Figure 2A).

On top of the interface, there is a toolbar. These additional tools allow you to have a better phenotyping experience:

- *Zoom In/Out:* you can zoom in and out on the display area (Figure 2B);
- *Pan:* after you zoom in, you may move the image at specific locations;
- *Data Cursor:* you can read a value at a specific pixel in the image.
- $! \rightarrow$ Many of these *Interface* arrangement are consistent across the other modules in Phenotiki.
- **Import Images:** In the bottom-left part of the window, press the button *Import* and it will be asked which import method you want to do. You have three different import options to add your images to the Pot Tray Analysis:
 - i. Images acquired from Phenotiki device;
 - ii. CSV File;
 - iii. Folder.

Images from device: The Phenotiki software integrates seamlessly with the device and can easily import images acquired using the Phenotiki device. When this import option is chosen, simply navigate to the folder containing the time-lapse images acquired with the Phenotiki device and press *Open*.

! → Make sure that image files conform to the following filename format: IMG_YYYY-MM-DD_HH-NN.png (cf. Table 2 for further details), which describes the acquisition date and time of that picture. This information is used within the software to automatically sort the time-lapse images according to their chronological order, as well as to extract time-dependent phenotyping traits (further details in Section 2.4).

CSV File: To import time-lapse images acquired with a device different than Phenotiki and you want to be sure that the temporal sequence is maintained, you need to build an index file first. This operation has to be done externally to the Phenotiki software.

For example, follow these steps:

- 1. open a spreadsheet software (e.g., Microsoft Excel);
- 2. fill the first column (A) with the filenames of your images;

Table 1. Frant traits extracted by the Fot Tray Analysis module.				
Trait	Unit	Description		
Projected Leaf Area	cm^2	Plant area calculated on the		
(PLA)		number of visible pixels		
Diameter	cm	Longest distance between two		
		points in the boundary		
Perimeter	cm	Length of the boundary pixels		
Compactness				
(or solidity)	-	Ranges in [0,1] indicating the		
		solidity of the plant (e.g., absence		
		of holes in the plant surface)		
Stockiness		-		
(or form factor)	_	Ranges in [0,1] indicating the		
		circularity of the plant		
Relative Growth Rate	$\%h^{-1}$	Measure of the amount of growth		
(RGR)		between to consecutive time		
		instants		
Color	0	Average value of the Hue (H)		
		component in the HSV color space		

Table 1: Plant traits extracted by the Pot Tray Analysis module



Figure 3: Expected output of the Pot Tray Analysis module.

- 3. fill the second column (B) with the timestamps referring to the time your images were taken (date format YYYY-MM-DD HH:NN, please refer Table 2 for further details);
- 4. enter as many rows as the images in the dataset (expected result is shown in Figure 4) and finally save your file in Comma-Separated Value format (CSV) *inside* the same folder where your images are located.

Then, in the Phenotiki software:

- 1. open the Pot Tray Analysis module;
- 2. press the "Load Dataset" button;
- 3. a new window for file selection will appear. In the bottom part, select as file type "*Comma-separated Values files (.csv)*";
- 4. navigate to the folder containing your images and the created CSV file;
- 5. choose the CSV file;
- 6. Press "Open".

Folder: With this option, you can import all the images contained inside a folder. Once this option is selected, you simply need to navigate to the desired folder and press "*Open*".

Table 2: List of abbreviation used for the image filename format.

Abbreviation	Meaning	Example
YYYY	Year (4 digits)	2016
MM	Month (2 digits)	01
DD	Day (2 digits)	14
HH	Hour (2 digits, 24hr format)	19
NN	Minutes (2 digits)	00

- $! \rightarrow$ Importing images with this method does not ensure temporal consistency, as the images will be added as they appear in the filesystem. We encourage to add them using CSV file instead.
- **Placing plant centers:** The segmentation algorithm implemented on Phenotiki relies on the prior knowledge of the position of plants [4]. In order to mark the position of each plant:
 - 1. Hold on the keyboard either:
 - On Mac OS X: 🕱
 - On other systems: ctrl
 - 2. as you hold the aforementioned key, move the mouse cursor on the first plant;
 - 3. press the left-button of your mouse to mark the plant;
 - 4. repeat these steps for all plants;
 - 5. when you have finished, release the key you have been holding.
 - $! \rightarrow$ For best results, place a plant's mark approximately on its centre and make sure that the mark is in the plant region.

As you add plant centers, the list in the right-hand side of the interface is filled with the coordinates of the clicked points. If plants belong to different group/genotype/treatment, you may specify it in the second column of the list of plants.

 $! \rightarrow$ Defining group/genotype/treatment will affect the leaf counting. Please read Section 2.3 for further details.

Parameter settings: The module exposes several parameters of the plant segmentation algorithm, that can be adjusted in case default values do not produce acceptable results. To set the parameters, you need to press the button "Settings"

	A	В	С
1	IMG_2013-09-28_08-00.png	2013-09-28 08:00	
2	IMG_2013-09-28_19-40.png	2013-09-28 19:40	
3	IMG_2013-09-29_08-00.png	2013-09-29 08:00	
4	IMG_2013-09-29_19-40.png	2013-09-29 19:40	
5	IMG_2013-09-30_08-00.png	2013-09-30 08:00	
6			
7			
8			
9			
10			
11			
12			
13			

Figure 4: Example of CVS for non-phenotiki dataset

Options	
Pixel size:	0.032258 cm
Number of plants:	24
Noise removal thr.:	60
Sigma:	1
Appearance Model:	
Lambda:	0.5 ∢ ▶
Features:	🗌 R 🔄 L 🔄 Gray 📄 LoG
	G 🗹 a* 🗹 TFB 📃 DoG
	B 🗹 b* 📄 Median 📄 Wiener
ĺ	Advanced Options Save

Figure 5: Window showing the set of parameter of the segmentation tool

and the window shown in Figure 5 will appear. Table 3 shows a brief list of the most important parameters and their meaning. Detailed explanation of the plant segmentation approach and all parameters can be found in [4].

	1 0	
Parameter	Description	Default value
Pixel size	Size of pixel in centimeters	Read below
Noise threshold	Remove objects having area smaller	60-80
	than a threshold (pixels)	
Sigma	Larger Sigma produces more	1
	regular boundaries, but thin petioles	
	could be missed	
Appearance Model	Use a plant appearance model	false
	to adapt to appearance changes	
Lambda	Influence of the appearance model	0.5
	(if used) compared to image	
	features (value in [0,1])	
Features	List of image features used	TFB, a*, b*
	during the segmentation	

Table 3: Parameters used in the plant segmentation tool.

- $! \rightarrow$ The *Number of Plants* is set automatically when providing plant centres.
- ! → Setting these parameters may require some expertise. Please go to Section 3.3 to learn more how to face frequently occurring issues.

Extract the masks: Once the centers of plants are placed and the parameters are set, you are ready to extract the segmentation masks, by pressing the button "*Extract Mask*". This operations will be applied to all the images that are ticked in

	the <i>File List</i> . On the right-hand side of the window, the list of plants with their location (in pixel) is shown. From here you can:
	• assign genotype/group name: by clicking on the corresponding cell, you can add/modify the name of the group of that plant. This operation is automatically propagated to all the images in the stack, assigning the same group to the same subject. Afterwards, you may edit group assignments in the time-lapse sequence.
	• view data : by pressing the right button of your mouse on the row of a plant a contextual menu will appear. There you can select <i>Properties</i> to get access to all the details of that particular subject. At this point, leaf counting is not provided. This will be discussed in Section 2.3.
Get the results:	Press on the button " <i>Get Traits</i> " to compute the phenotyping traints in Table 1. Once plant segmentation is done, results are displayed as in Figure 3.
! ightarrow	Save your progress! So far your data are not stored on your hard drive. By pressing the button " <i>Save</i> ", you can save your data. We recommend to do it at each step of your analysis to avoid any data loss. The database of your data is a <i>.mat</i> file. To learn how to extract plant traits, go to Section 2.4.
Modality selection:	On the right-hand side of the interface, you can find a pop-up menu allowing to switch between different visualization modalities. For this module, you can visualise:
	1. <i>raw image</i> : the original RGB images;
	2. detected plants: colour-coded segmentation masks (cf. Figure 3);
	3. <i>contour:</i> RGB images with contour overlaid on the plants. This modality is helpful to assess visually the goodness of the segmentation;
	4. <i>FG mask:</i> segmentation mask where foreground (plants) is white and background is black.
Load your data:	If you wish to load data from a previous analysis, press the " <i>Load Dataset</i> " button and navigate to the folder containing the database file with <i>.mat</i> extension.
2.2 Leaf Annotation	
	This module is largely based (and improves upon) our stand-alone Leaf Annotation tool, for which a detailed manual is available on the Phe- notiki website (http://phenotiki.com/download.html). The leaf segmentation engine is based on the random walker algorithm [5] and the overall approach is described in [3].
	The user annotate the image by drawing scribbles on the leaves (at least a scribble per leaf, for all the leaves in the plant). The suggested workflow for this module is the following: <i>load dataset, leaf annotation</i> , β <i>selection, segmentation</i> .
Load Dataset:	On the left-hand side of the interface, press on the " <i>Load Dataset</i> " button. Once the dialog window is open, navigate to the folder containing the <i>.mat</i> you prepared beforehand (e.g., from the Pot Tray Analysis). The interface allows you to browse through tray and plant images.
Leaf annotation:	Once you have selected the plant to annotate, you will need to draw a scrib-

ble for each leaf. The software will use these scribbles to delineate each individual leaf and generate a leaf segmentation of the plant. Scribbles can be of several types:

- *Dot:* a single pixel is marked to belong to a certain leaf;
- *Line:* all the pixels laying on a line segment are marked to belong to a leaf;
- *Freehand:* a freehand line is marked to belong to a specific leaf.
- ! → If you feel more comfortable with using other third-party image manipulation software (e.g., Photoshop, GIMP), you can also use them to create a binary image for an annotation. You can import such external scribbles by pressing the "*Load from file...*" button.

As you add scribbles, on the right-hand side a list of annotations will be populated. They are organised hierarchically:

- i *label:* a label is a container of scribbles identifying a single leaf;
- ii *scribble:* one or more scribbles are contained in a label, and each will mark pixels that belong to the same leaf.

In many cases, a single dot scribble is enough to obtain a suitable leaf segmentation. However, mature plants with heavily overlapping leaves may require more detailed annotations to produce an accurate leaf segmentation.

From the list of labels and leaves, you can access a contextual menu (using the right-button of your mouse) allowing you to:

- delete labels;
- delete scribbles;
- change label to a scribble;
- copy labels applied to the same plant at the previous time instant.
- $! \rightarrow$ Each label is assigned a numerical id. When you delete a label, the numerical id of the others will not change. For example, if you have *Label 1*, *Label 2*, and *Label 3*, and you delete *Label 2*, than *Label 3* will not change. This is to maintain label consistency across time.
- β selection: This is the only parameter required for this module. It controls the sensitivity of the segmentation algorithm to variations in colour intensity between neighboring image pixels. Higher values of β correspond to lower sensitivity to sudden changes in color intensity. Example values are: 15, 30, 80. Experiment with different values to find the one that produces best results on your dataset.
- **Leaf segmentation:** Once labels have been placed, press the "*Segment*" button in the lower part of the interface. The time to required by the software to generate the segmentation is usually in the order of few seconds and depends on image resolution. When it is done, the software will show you the outcome.
 - **Refine result:** If you are not satisfied with the outcome, you may add further scribbles to obtain refined results. For example, you can change a dot annotation into a line, or even replace it with a more detailed freehand scribble. In some cases, it could be useful to add multiple scribbles to the same leaf label. To add multiple annotations to the same label, go to the list of labels

on the right-hand side of the interface and press on the desired label (e. g., *Label 1*). At this point, choose one of the tools (e. g., line) and annotate as before. This new scribble will appear with the same colour of the other scribbles within the same label. If you want add more labels, press on '*<new label>'* from the same list.

! → Note that the number of labels can be used by the *Leaf Counting* module, which will be discussed in the next section. Therefore, if you are more interested in the number of leaves than in obtaining a leaf segmentation, you can use this module to produce "training data" for the learning-based counting algorithm. In this case, using the dot annotation tool is an efficient method to rapidly annotate the number of leaves for several plants.

2.3 Leaf Counting

The Leaf Counting module is able to count the number of leaves in a plant image. It is based on a *machine learning* algorithm, so it can adapt to different types of plants and varying imaging conditions. It can do so by relying on a sufficient amount of labelled data. Therefore, the user has to provide leaf count for several plants, such that the software can learn a model that will be applied to the rest of the dataset.

The suggested workflow for this module is the following: *data loading*, *model training*, *testing*.

- **Interface:** In Figure 6 the Leaf Counting module interface is shown. It has several panels and they will be referenced next in the manual. A brief explanation is given as follows:
 - A. **Plant list:** here you can load data that you previously processed with the other modules. The interface lists all plants across all images and provides a summary showing: total number of plants (i.e. number of subjects multiplied by number of time-lapse pictures), training set size (i.e. number of labelled data), and testing set size (i.e. number of unlabelled data).
 - B. Learning setup: here you set which plants will belong to the training and testing set, respectively. Although only labelled plants will be included in the training set, you have the possibility to filter them by group (e.g., to test performance on a specific group of plants). You can also assign a name to a model, such that you can train several models (e.g., corresponding to different training sets or parameters) according to specific needs.
 - C. **Results:** every time a model is learnt, a new entry on the results list will appear. From there, you may check the performance of the learning algorithm, as well as repeat the same experiment with different parameters. When a leaf count is estimated, these values are not automatically applied to the dataset. In order to apply leaf counting estimates to the dataset press on the *Assign Count* button.

Each learned model will appear separately in the results list.

In addition, by pressing on *Show Advanced Options...*, the following panels are also displayed (not showed in Figure 6):

• **Model parameters:** the set of parameters for the learning algorithm can be tuned from this panel.

Sequence: 28-Sen-2 *				
Sequence. Eo Sep Ent	Trainin	g samples: All labelled i	mages	
Train Test # Group Leaves	Testin	g samples: All unlabelle	d images	;
	Name of the test: Noname			
8 ein2.1 9 ctr1		Dataset loade	Dataset loaded successfully	
	Start Training		Show Advanced Options	
12 ctr1	Start Testing		Save	
✓ 13 pgm ✓ 14 adh1	Fresults			
🗌 🗹 15 ein2.1	Noname			
16 ctr1				
17 adn1				
19 ctr1				
20 Col-0				
				_
	Metric	Value		
Total number of plants: 120	1 CountDiff	+0.00000 (0.00000)		
Total number of plants: 120	a Absolute CountDiff	10.00000 (0.00000)		
Total number of plants: 120 Training set size: 1	2 Absolute CountDiff 2 Percent Agreement	+0.00000 (0.00000) 100.00%		
Total number of plants: 120 Training set size: 1 Testing set size: 119	2 Absolute CountDiff 3 Percent Agreement 4 Mean Squared Error	+0.00000 (0.00000) 100.00% 0.00000		
Total number of plants: 120 Training set size: 1 Testing set size: 119	2 Absolute CountDiff 3 Percent Agreement 4 Mean Squared Error 5 R-squared	+0.00000 (0.00000) 100.00% 0.00000 NaN		
Total number of plants: 120 Training set size: 1 Testing set size: 119 Load Dataset	Absolute CountDiff Apsolute CountDiff A Peccent Agreement Mean Squared Error S R-squared	+0.00000 (0.00000) 100.00% 0.00000 NaN	Logpolar	•

Figure 6: Leaf counting module. Panels are explained in text.

- **Cross-validation:** in order to check whether your model is well-trained, you may run a cross-validation. It will split the training set into K folds, out of which K 1 are used for training, and the remaining one is left for testing. K rounds of validation will be performed and average results will be reported.
- **Display:** training and test counts are displayed. Moreover, performance plot for training error may be shown.
- **Load data:** Press on the "*Load Dataset*" button to load data you previously processed with the other modules. The number of count used for the training set can be provided using the following methods:
 - Leaf Annotation Tool: using the leaf annotation tool, you will also provide the number of leaves in all the annotated plants. This information will be used by the Leaf Counting module as training data.
 - **Manual prompt:** the list of subjects in the Plant List panel (cf. Figure 6) has 5 columns. The last one contains the number of leaves that the software already found in the database (either by previous counts, or leaf annotation tool). By double-clicking on a cell you can input a value.
 - External file: you may load the number of leaves using an external CSV file. It is a similar import procedure as discussed in Section 2.1. In an external spreadsheet software (e.g., Excel), create a column for each subject in the image and record the known leaf count for each time instant. Figure 7 shown visually how to format the CSV file. Please remember to save the file in *.csv* format.
- **Model training:** Once a certain amount of training data is available (i. e. plants for which the number of leaves is known), you can train the regression model. Firstly, select the training samples (default is all labelled plants, but for example



Figure 7: Loading known leaf counting from a CSV file. Using a spreadsheet program (e.g., Excel), you may create your counting file to provide labelled data to the leaf counting module. Each column corresponds to a plant subject, whereas the rows correspond to time vertically you record the number of leaves across time. Plants and columns have to match with respect to Plant ID added in the Pot Tray Analysis module (bold yellow number). Unknown entries have to be marked with **zeroes** '0'. Remember to save the final file in *.csv* format.

you may train only on specific genotypes), then press on **Start Training**. Once this operation is completed results will be displayed (Figure 6C). Please refer to [6] for further explanation about the performance measures.

- $! \rightarrow$ If the *Advanced Options* CategoricalField is defined with Group (which is done by default), you must add a plant for each group/genotype/treatment you have defined. If no group/genotype/treatment was set with the Pot Tray Analysis, you can ignore this alert.
- $! \rightarrow$ The time required for this operation depends on the size of the training set. It may require several minutes to be completed. In general, more training data lead to better performance.
- $! \rightarrow$ You may train several models according to different training samples or different parameters. We advice you to keep track of your changes by providing a unique name to each test. Tests with same name will be overwritten.
- **Testing:** To count leaves of plants with unknown counting, select a trained model from panel *C*, select the testing set form panel *B*, and finally press **Start Testing** button.
 - $! \rightarrow$ At this stage, the leaf counting is **not yet applied** to the subjects, because you may want to do several tests before applying estimated counts. Once you are satisfied with the count, press **Assign Count** button in panel *C*.
- **Parameters:** The leaf counting model may require to the user to set some parameters. A full list of parameters with description is provided in Table 4. Interested readers may refer to [6] for further details. To test the performance of a set of parameters, we suggest to run *K*-fold cross validation. You may run a cross-validation pressing on the *Show Advanced Options...* button.

Parameter	Explanation	Suggested value	
Patch size	Size of image patches from images that are extracted	1/4 of the average-sized leaf	
Ratio Curve	Size of the sliding window in the	As wide as an	
Window Width	log-polar domain to detect patch candidates	average-sized leaf	
LogPolar Pooling Bands	Number of areas where features are pooled together in the log-polar space	5	
Cartesian Pooling Bands	Number of areas where features are pooled together in the Cartesian space	3	
Pooling Callback	Function used for the pooling	max (max pooling)	
LogPolar	Rescaling value in order to	0 (dynamic scaling)	
Normalization	normalize the height of log-polar images		
Cartesian	Rescaling value in order to	100 (100x100 pixels)	
Normalization	normalize the size of Cartesian images		
Cartesian Features	Boolean value: 1 means that Cartesian features are used, 0 otherwise	1	
LogPolar Features	Boolean value: 1 means that log-polar features are used, 0 otherwise	1	
Dictionary Size	Number of clusters for the k-means, which builds the patches dictionary	50	
Pooling Shifts	Dataset augmentation via pooling shifts	4	
SVM Matrix	Normalization function for the SVM data	zscore	
Normalization			
SVR_OPT.C	SVR Optimized parameters – Loss term	1	
SVR_OPT.gamma	SVR Optimized parameters – RBF gamma	0.003	
SVR_OPT.epsilon	SVR Optimized parameters – Allowed regression error	1	
RF.NumTree	Random Forest – Number of trees	100	
RF.LeafMinSize	Random Forest – Minimum number of data points per leaf	5	
Time Features	Adds a variable describing day a plant was acquired (date and time must be set)	Left blank	
Categorical Field	Property of the plants used as categorical variable	Group	

Table 4: Description of the parameters of the leaf counting algorithm.

2.4 Data Extraction

This module allows you to extract plant traits from the image data analysed using all other modules. You may extract the following data:

- Plots: you can plot phenotyping traits across time and save the plot;
- **Images:** you can save images from the dataset (e.g., segmentation masks, annotated leaves);
- Raw data: data extracted from plants are formatted into a CSV file

6000 -		adh etri
4000		
		pgn
2000 -		
0000 -		
8000 -		
6000 -		
4000 -		
2000		
0 28/Sep	29/Sep	30/Sep
	2000	2000 000 000 000 000 000 000 000

Figure 8: Data extraction module.

for subsequent analysis in external tools.

2.4.1 Data plot

Once you load your dataset, you can plot data of the phenotyping traits extracted from your plant images. The list of data available for plotting is shown in Table 1. An example of the interface is depicted in Figure 8.

When you plot the data for a specific traits, you may set the span of time. You also have three different modalities of plotting: (i) all, (ii) specific group, and (iii) specific subject. By selecting *all*, the software will plot the mean value as a line of the trait for each group (please go to Section 2.1 to learn how to associate plants to groups); the line will appear in the middle of a wider band, indicating the lowest and higher bound for that trait. Instead, by selecting *specific group*, the software will plot all the subjects individually that belong to a specific group/genotype. You may change the group from the drop-down menu *specify*. By selecting *specific subject*, the software will plot the selected trait of the specified subject. Plots can be exported by pressing **Save Plot As**. They will be saved as images in PNG format.

 $! \rightarrow$ To convert the measures in centimeters, you need to tick the box *Converts* values in cm (cf. Figure 8). You must specify the value *Pixel size* from the Settings in the Pot Tray Analysis tool (cf. Figure 5).

2.4.2 Data and images extraction

By pressing on *Extract Data*, you launch the export wizard to extract data or images.

Extracting data. To extract data, select the item *Data* and then press **Next**, then select which traits you want to extract. You can choose multiple items (press **Select all** if you want all of them) and then press **Next**. The software will ask you to choose the folder where data will be saved.

 $! \rightarrow$ Each trait will generate a different CSV file.

Extracting Images. To extract images, choose the item *Images* and then press **Next**. You can extract either *Individual Plants* or *Entire Trays*. Once you have done your choice, press **Next** and select the folder where you want to save the images.

3 Troubleshooting

3.1 General

Phenotiki is telling me that *libsvm* was successfully installed, but I get some errors in the Matlab console. How can I fix it? The installation script within *libsvm* has an error handle functionality that cannot be overridden by Phenotiki. Often, this problem appears when you do not have a suitable C compiler installed and/or configured in your system. We recommend to read the follow instruction at https://uk.mathworks. com/support/sysreq/files/SystemRequirements-Release2015b_ SupportedCompilers.pdf

How can I move the data to/from another workstation? If you wish to share data with colleagues or move the data to another workstation, you simply need to copy the plant's image files and the database file with *.mat* extension that was generated by the modules. When you load your data on the new workstation, the software will ask you to provide the path where images are located.

 $! \rightarrow$ This procedure is applicable also in case you move the folder of the images somewhere else on the same computer.

Some computations take too long. Is there any way to stop them? Unfortunately it is not possible to abort a computation. Since our software is based on MATLAB and it is single-threaded, we could not add any mechanism to cancel long computations. However, if you are running Phenotiki from the Source Code, you may press CTRL+Z, as you would do for any computation in MATLAB you want to abort.

- 3.2 Operating system-dependent issues
- 3.2.1 Windows

Problems with *mclmcr.dll*: If you have errors saying *mclmcr.dll not found*, it could be either of the following: (i) MCR not installed, (ii) MCR installation corrupted, (iii) Older version of the library installed in a different directory. Specifically, we suggest:

- i. Download the Installer of Phenotiki and make sure you are installing MCR. Alternatively, if you have the executable, follow the instruction how to install MCR in Section 1.2.
- ii. Different solutions can be found on the Internet (e.g. https://uk. mathworks.com/matlabcentral/answers/101146-whydo-i-receive-an-error-message-of-mclmcr-dll-notfound-while-deploying-my-application-on-a-targ. Search for MCR on your system and delete all the files related to it. Then, go to the point (i).
- iii. Search on your system for the *mclmcr.dll*. It can happen it is located within the System32 directory of Windows. Locate and remove it. Then, go to the step (ii).

Problem extracting zip files: If you are experiencing the error *An error occurred while extracting files* and you have the zip file containing Phenotiki in your Desktop folder, try to move the zip file in any other folder different than the Desktop and try again. We have experienced this issue and we do not know why it is happening.

Using virtual machines: We experienced issues using MATLAB/Phenotiki on Parallels (for Mac). We do not recommend using Parallels with Phenotiki on a virtual Linux Operating System.

3.3 Pot Tray Analysis

Cannot find suitable parameters. What can I do? Finding good parameters for the plant segmentation may not be immediate for challenging datasets. If you have fair results and you want to improve them, then we suggest to select a few images to trial different configurations. In the unluckily case that segmentation results are not accurate enough, we suggest to run a *grid search*. With this method, the software will automatically test a range of parameters and find optimal ones according to some supervision.

First, segment manually a few of plant images using an external image editing software (e. g., GIMP, Photoshop). Create a black-and-white (also known as binary) image with just 0's and 1's, where 0 (black) means background and 1 (white) means plant. When you are done with this, open the Phenotiki Analysis Software and go to Pot Tray Analysis tool and load your dataset. Then, from the list of images, select one plant image for which you have manually created the foreground/background segmentation and select it by with the right-button of your mouse. A contextual menu will appear and press *Set ground-truth*. When the file selection window appears, select the corresponding mask for this image and repeat this operation for all the plant images for which you have obtained the segmentation mask.

Then press on the **Setting** button in the bottom part of the Pot Tray Analysis Tool and press on **Advanced Options** (cf. Figure 5). The full list of parameters will appear and you may choose the ones you want to include in the grid search, by ticking the check-boxes in the list. When you choose a parameter, you may define a range of values that the software will test against the ground-truth you have provided before. Ranges can be set in the following way.

- numeric values: MATLAB syntax is accepted. For example:
 - explicit list of values to test: [0, 10, 20, 30, 40, 50];
 - all the values between 0 and 50 with a step of 10: 0 : 10 : 50;
- **text values:** cell array MATLAB syntax is accepted. Some examples are:
 - {Feature1, Feature2}: the software will use Feature1 and Feature2 separately;
 - {Feature1,Feature2}, {Feature1,Feature3}: the software will use Feature1 and Feature2 at the same time, and then Feature1 and Feature3 at the next test;

- **Boolean values:** for Boolean values you do not need to specify range of values, as they can be either 0 or 1.
- $! \rightarrow$ The more images have the ground-truth, the better it is in terms of accuracy. Notice though that the grid search will take longer to complete.

3.4 Leaf Annotation

A leaf is split in half by another leaf. What should I do? In this case, we suggest to add (at least) two annotations for the same leaf. Once you annotate one part of a leaf and a new label is created, select the newly created label from the list to keep adding annotations. You will have two or more annotations with the same label. To annotate a new leaf, select <*new label>* from the list of label.

I would like to keep consistent labels across time. If you want to assign the same label to the same leaf across time, we recommend to start the annotation from the first time instant (or the last). Then, move to the next time instant and from the list of labels click with the right-button of your mouse and select *copy from the previous plant* (or *copy from the next plant*). In this way, annotations are copied with consistent numeric identifiers. Then, edit the copied labels to adjust to growth and movement of the plant and add new labels if necessary.

Part of a leaf is assigned to another leaf. Use more refined annotation (such as lines or freehand scribbles) to better delineate the image region of each leaf.

3.5 Leaf Counting

What are the allowed values for the paramters? Whilst integers paramters may accept any value, we suggest you to vary them by small amounts. For the other parameters, please refer to Table 5.

Table 5:	Descriptio	on of the p	arameters	of the leaf	counting	algorithm

Parameter	Allowed Values
Pooling Callback	<i>max</i> or <i>sum</i>
SVM Matrix Normalization	<i>zscore, range0:1, range-1:1, norm2, none</i>
Time Features	Empty, <i>delta, unix</i>

Copyright and License

3.6 Phenotiki License

Copyright (c) 2016 Massimo Minervini, Mario Valerio Giuffrida, Sotirios A. Tsaftaris All rights reserved.

This license granted to You for the Phenotiki Analysis software (hereafter termed Licensed Application) by the authors is limited to a nontransferable license to use said software. You may not distribute or make the Licensed Application available outside your lab or institute. You may not rent, lease, lend, sell, redistribute or sublicense the Licensed Application. You may not copy, decompile, reverse engineer, disassemble, attempt to derive the source code of, modify, or create derivative works of

	the Licensed Application, any updates, or any part thereof (except as and only to the extent any foregoing restriction is prohibited by applicable law or to the extent as may be permitted by the licensing terms governing use of any open sourced components included with the Licensed Applica- tion). Any attempt to do so is a violation of the rights of the Application Provider and its licensors. If You breach this restriction, You may be subject to prosecution and damages. The terms of the license will govern any upgrades provided by Application Provider that replace and/or sup- plement the original Product, unless such upgrade is accompanied by a separate license in which case the terms of that license will govern. There is no warranty for the program, to the extent permitted by ap-
	plicable law. The program "as is" without warranty of any kind, either expressed or implied, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose. the entire risk as to the quality and performance of the program is with you. Should the program prove defective, you assume the cost of all necessary servicing, repair or correction.
vl_feat license	
	Copyright (C) 2007-11, Andrea Vedaldi and Brian Fulkerson Copyright (C) 2012-13, The VLFeat Team All rights reserved.
	Redistribution and use in source and binary forms, with or without mod- ification, are permitted provided that the following conditions are met:
	• Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
	• Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.
	THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PAR- TICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT HOLDER OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CON- SEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CON- TRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.
libsvm license	
	Copyright (c) 2000-2014 Chih-Chung Chang and Chih-Jen Lin All rights reserved.
	Redistribution and use in source and binary forms, with or without mod- ification, are permitted provided that the following conditions are met:

	• Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
	• Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.
	• Neither name of copyright holders nor the names of its contributors may be used to endorse or promote products derived from this software without specific prior written permission.
	THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WAR- RANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WAR- RANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICU- LAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE RE- GENTS OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAM- AGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUB- STITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.
findjobj license	
	Copyright (c) 2015, Yair Altman All rights reserved.
	Redistribution and use in source and binary forms, with or without mod- ification, are permitted provided that the following conditions are met:
	• Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
	 Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribu- tion
	THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PAR- TICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT OWNER OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CON- SEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CON- TRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

Boundedline license

The MIT License (MIT)

Copyright (c) 2015 Kelly Kearney

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

References

- [1] A. Vedaldi and B. Fulkerson, "VLFeat: An open and portable library of computer vision algorithms," 2008, http://www.vlfeat.org/.
- [2] C.-C. Chang and C.-J. Lin, "LIBSVM: A library for support vector machines," ACM Transactions on Intelligent Systems and Technology, vol. 2, pp. 27:1–27:27, 2011, software available at http://www.csie. ntu.edu.tw/~cjlin/libsvm.
- [3] M. Minervini, M. V. Giuffrida, and S. A. Tsaftaris, "An interactive tool for semi-automated leaf annotation," in *Proceedings of the Computer Vision Problems in Plant Phenotyping (CVPPP) Workshop.* Swansea, UK: BMVA Press, Sep. 2015, pp. 6.1–6.13.
- [4] M. Minervini, M. M. Abdelsamea, and S. A. Tsaftaris, "Image-based plant phenotyping with incremental learning and active contours," *Ecological Informatics*, vol. 23, pp. 35–48, Sep. 2014.
- [5] L. Grady, "Random walks for image segmentation," *IEEE Transactions on Pattern Analysis and Machine Intelligence*, vol. 28, no. 11, pp. 1768–1783, Nov. 2006.
- [6] M. V. Giuffrida, M. Minervini, and S. A. Tsaftaris, "Learning to count leaves in rosette plants," in *Proceedings of the Computer Vision Problems in Plant Phenotyping (CVPPP) Workshop.* Swansea, UK: BMVA Press, Sep. 2015, pp. 1.1–1.13.